

Experiment AMe-1: Small Animal Respiratory Exchange Ratio (RER)

Equipment Required

PC or Mac Computer

iWorx IXTA, USB cable, and power supply

Small Animal Chamber

Gas sample tubing

GA-200 CO₂/O₂ gas analyzer with filter

A-CAL-150 Calibration Kit

Barometer and thermometer

Mouse or Rat

Snake, Frog or Medium sized Lizard

Balance or scale

Setup the Metabolic Cart

1. Locate the small animal chamber, the GA-200 gas analyzer, the gas analyzer power supply, two sensor output cables, a gas inlet filter, and Nafion sampling tubing.
2. Position the gas analyzer on the desktop, so that the analyzer can be connected to the data recording unit and the small animal chamber at the same time.
3. Place the clear end of the Nafion sampling tube on the gas sampling port near the outlet of the small animal chamber.
4. Place the filter on the inlet port in the lower right front corner of the gas analyzer. Attach the braided end of the Nafion sampling tube to the filter.




Figure AMe-1-S1: The small animal chamber with an endothermic animal inside.

5. Connect the outputs of the oxygen and carbon dioxide sensors, which are located on the rear panel of the gas analyzer, to the BNC inputs of the recorder.



Figure AMe-1-S2: The front panel of the IXTA.



Figure AMe-1-S3: The GA-200 gas analyzer showing the proper connections of the Nafion tubing to the front filter to measure O_2 and CO_2 exchange.

Set up the GA-200:

- Connect the BNC output of the carbon dioxide sensor, labeled CO2 SENSOR, 0.8-4V=0-10%, to Channel A3 using a male BNC-BNC cable.
- Connect the BNC output of the oxygen sensor, labeled Oxygen Sensor and 0-1V=0-100% to Channel A4 using a male BNC-BNC cable.



Figure AMe-1-S4: The rear panel of the GA-200B gas analyzer showing the outputs of the carbon dioxide sensor, on the left, and the oxygen sensor, on the right.

7. Plug the GA-200 into an outlet.
8. Use the power switch to turn on the gas analyzer. As the unit powers up, the display on the front of the unit will light up.

Calibrating the Gas Analyzer

The gas analyzer is calibrated and checked at the factory and has the ability to measure and record both oxygen and carbon dioxide concentrations over the ranges normally recorded from human and animal subjects.

However, prior to use, users must perform a two-point calibration of the oxygen and carbon dioxide sensors.

The recommended gas mixtures include:

- 12% O₂ and 5% CO₂ with the balance of the mixture being N₂; or,
- 16% O₂ and 4% CO₂ with the balance of the mixture being N₂.

Room air can provide both the high concentration of O₂ at 20.90%, and the low concentration of CO₂, 0.04%.

Note: Warm up the gas analyzer for at least 30 minutes. The input Air Filter must be connected to avoid damaging the sensors.

Record the Voltage Outputs of the Gas Sensors

1. Turn on the gas analyzer for at least 30 minutes before performing a calibration.
2. Prepare the equipment that will deliver any gas samples, other than room air, to the GA-200:
 - Clamp and secure any gas cylinders that will be used to provide gas samples near the gas analyzer.
 - Attach the regulator to the gas cylinder safely.

- Attach a Luer-Lock connector to the outlet of the regulator that will allow the Calibration Kit for the gas analyzer to be connected to the regulator of the gas cylinder.
 - Open the primary and secondary valves of the regulator for a few seconds to purge the air from the regulator.
 - Close the secondary valve on the regulator to stop the flow of gas from the regulator. You will need the cylinder for the second sample of gas.
3. Attach a filter to the inlet port on the gas analyzer. Attach the braided end of the sampling tubing to the inlet of the filter.
 4. Attach the clear tubing from the outlet port on the gas analyzer to the small animal chamber.
 5. Measure the voltage outputs of the oxygen and carbon dioxide sensors when measuring a sample of room air.
 6. Place the gas sampling tubing away from the users to prevent the sampling of exhaled air. Allow room air to be pumped through the gas analyzer for 10 seconds before recording the outputs of the sensors.
 - Type **Room Air** in the Mark box to the right of the Mark button.
 - Click on the Record button. The recording should scroll across the screen.
 - While recording, press the Mark button to mark the recording with information about the room air gas sample.
 - Record the outputs of the two gas sensors for about ten seconds. The recording which should be like the first segment of data below.
 - Continue to record while moving to the next series of steps.
 6. Measure the voltage outputs of the oxygen and carbon dioxide sensors when measuring a second sample of a gas mixture containing known concentrations of oxygen and carbon dioxide.
 - Connect the A-CAL-150 Calibration Kit to the Luer-Lock connector on the output of the regulator.
 - Connect the A-CAL-150 or 200 Calibration Kit to the inlet filter port on the gas analyzer using the nafion tubing. The gas analyzer will pull the air in from the calibration kit.



Figure AMe-1-S5: Calibration Kit A-CAL-150

- Open valve on the regulator of the cylinder providing the gas sample.
 - Type **Gas Sample** in the Mark box.
 - While continuing to record with the sample gas flowing into the gas analyzer, press the Mark button to mark the recording with information about the gas sample.
7. Once the recordings of the gas concentrations reach a steady level, record for another ten seconds.
 8. This may take a minute or so.
 9. Click the Stop button.
 10. Select Save As in the File menu, type a name for the file. Click on the Save button to save the data file.

Convert the Units on Gas Concentration Channels

1. Use the Display Time icons to adjust the Display Time of the Main window to show the complete calibration data on the Main window at the same time. The required data can also be selected by:
 - Placing the cursors on either side of data required.
 - Clicking the Zoom between Cursors button on the LabScribe toolbar to expand the entire segment of data to the width of the Main window.




Figure AMe-1-S6: The voltage outputs of the two sensors in the gas analyzer, carbon dioxide on the top and oxygen on the bottom. Other recording windows have been minimized to show detail.

- Click the 2-Cursor icon so that two blue cursors appear on the Main window. Place one cursor on the section of data recorded when gas analyzer was collecting a sample of room air and the second cursor on the section of data recorded when the second sample was collected.




Figure AMe-1-S7: The LabScribe toolbar.

- Convert the voltages at the positions of the cursors to concentrations using the Advanced Units Conversion dialogue window.
- To convert the voltages on the Expired CO₂ Concentration (%) channel, click on the arrow to the left of the channel title to open the channel menu. Select Units from the channel menu, and select Advanced from the Units submenu.




Figure AMe-1-S8: The Advanced Units Conversion dialogue window with the voltages between the cursors set to equal the concentrations used in calibration.

6. On the Units Conversion window:

- Make sure Apply units to the next recorded block and Apply units to all blocks are selected in the menu under the displayed graph on the left side of the window by putting a check mark in the boxes next to each statement.
- Click on and move the cursors so that they are in position such that:
 - the first 2 cursors are on the area where room air values were recorded. Leave a space between the cursors so that you have an average value being calculated while room air was moving into the GA-200 gas analyzer.
 - the second 2 cursors are on the area where the gas sample values were recorded. Leave a space between the cursors so that you have an average value being calculated while the gas sample was moving into the GA-200 gas analyzer.
- Notice that the voltages from the positions between the cursors are automatically entered into the value equations. Enter the two concentrations of carbon dioxide measured from the two samples in the corresponding boxes on the right side of the conversion equations.
- Using room air - the concentration of CO₂ = 0.04
- The second gas concentration will be the one from the gas cylinder. Generally a 5% CO₂ concentration is recommended.
- Enter the name of the units, %, in box below the concentrations.
- Click the OK button in the lower right corner of the window.

7. Repeat Steps 4 and 5 on the Expired O₂ Concentration (%) channel.

- Room air = 20.9
- Second gas concentration will be the one from the gas cylinder. Generally a 12% O₂ concentration is recommended.

8. Click on the Save button.

Note: When using this LabScribe calibration protocol, the numbers in the software do not correlate with those on the front panel of the GA-200 gas analyzer.

Experiment AMe-1: Small Animal Respiratory Exchange Ratio (RER)

Exercise 1: Changes in CO₂ and O₂ in a closed chamber, and RER in an endothermic animal.

Aim: To determine mean RER of an endotherm at rest.

Approximate Time: 15-20 minutes or longer depending on the animal

Procedure

1. Place the animal in the chamber and close the lid securely.
2. Click on the Record button. Type **Mouse/Rat** in the Mark box to the right of the Mark button. Press the Mark button to mark the recording.
3. Click the AutoScale All button.
4. On the Expired CO₂ Concentration (%) channel, notice that the CO₂ concentration shows a continuous steady rise in level throughout the recording.

Note: During the first few minutes or so of the recording, the small animal chamber is filling with expired air from the organism inside. Since this is a closed system, the concentration of CO₂ continues to rise as the animal exhales into the chamber.

- The time that it takes the chamber to be filled with expired air and a steady increase in the level of carbon dioxide will depend on the volume and respiration rate of the organism and the volume of the small animal chamber. It will take longer to fill the chamber if the organism's respiration rate and tidal volume are low, or the animal is very small.

Warning: The CO₂ concentration will continue to rise throughout the experiment and levels can become toxic to the animal in the chamber. Remove the animal immediately if the CO₂ levels reach 3.5%.

5. On the Expired O₂ Concentration (%) channel, notice that the O₂ concentration will show a steady decrease as the animal is using the O₂ in the chamber. As pointed out in the previous step, the size of the chamber, the tidal volume, and respiration rate of the organism, determine the time it takes for the concentration of oxygen to decrease significantly.
6. On the RER channel, an equation programmed in the software determines RER based on the relative levels of CO₂ (which is increasing) and O₂ (which is decreasing).
7. Continue to record until the concentration of carbon dioxide in the chamber reaches 3%.

Note: Remember to remove the animal from the chamber immediately if the carbon dioxide concentration reaches 3.5%.

8. Once the appropriate data is recorded, click Stop. Your recording should be similar to the data displayed below.
9. Select Save As in the File menu, type a name for the file. Click on the Save button to save the data file.
10. Remove the animal from the chamber and place it back into its container.
11. Clean and dry the small animal chamber if necessary.

RER Channel Set Up

1. Display the complete data recording in the Main window. Use the Display Time icons to adjust the Display Time of the Main window to show the complete recording on the Main window.
2. Select and display a section of data that shows the change in CO₂ concentration at 1% as shown below.
3. This can be done by:
 - Clicking on and dragging the cursors to either side of the data, and looking at the V2-V1 value in the upper right corner of the Expired CO2 channel.
 - When V2-V1 value is equal to 1%, look at the T2-T1 value in the upper right to determine the time it took for the CO₂ concentration to change by 1%.
4. Click on RER Expired CO2 Concentration (%) on the RER channel. Choose Setup Function from the drop down list. This will open the RER Calculation Dialog window.
5. Setup the RER calculations:
 - Choose the appropriate RER Type by clicking on the down arrow and choosing Closed Small Animal Chamber.
 - Check the O2 and CO2 channel information so that they are being calculated from the correct channels.
 - Change the Time(s) to average to 30 seconds
 - Change the Delta Time(min) to be the T2-T1 value noted in Step 2. This is the time it took for the CO2 concentration in the small animal chamber to change by 1%.
 - Click OK.
6. The recording on the main window will now have a histogram on the RER channel.




Figure AMe-1-L1: The recording of the expired CO₂ over a period of 28 minutes showing a change in concentration by 1%.

Figure AMe-1-L2: RER Calculation Dialog window.

Data Analysis

1. Scroll to the beginning of the recording where the histogram is not a straight line. You should be able to see a step-like histogram.
2. Click and drag the cursors to either side of the step-like histogram.

- Click the Zoom between Cursors button on the toolbar to expand this section of data to fill the window.
- Click AutoScale on all channels. Your data should look like the data seen below.




Figure AMe-1-L3: The RER channel showing the histogram for the ratio of CO₂ to O₂ in the endothermic organism.

- Click on the Analysis window icon in the toolbar or select Analysis from the Windows menu to transfer the data displayed in the Main window to the Analysis window.
- Look at the Function Table that is above the uppermost channel displayed in the Analysis window. The functions, Mean, T2-T1, and V2-V1 should appear in this table. Values for these parameters on each channel are seen in the table across the top margin of each channel.
- Once the cursors are placed in the correct positions for determining the CO₂ and O₂ concentrations, the values for these parameters can be recorded in the on-line notebook of LabScribe by typing their names and values directly into the Journal.
- The functions in the channel pull-down menus of the Analysis window can also be used to enter the names and values of these parameters in the Journal. To use them:
 - Place the cursors at the locations used to measure the concentrations and RER.
 - Transfer the names of the mathematical functions used to determine these data to the Journal using the Add Title to Journal function in the Expired CO₂ Channel pull-down menu.
- Transfer the values for the data to the Journal using the Add All Data to Journal function in the Expired CO₂ Channel pull-down menu.

10. Repeat this procedure for the Expired O₂ Concentration (%) channel.
11. The values for the following parameters are determined when the cursors are positioned as directed:
 - Mean concentration of CO₂ in expired air, which is the value for Mean on the Expired CO₂ Concentration channel.
 - Mean concentration of O₂ in expired air, which is the value for Mean on the Expired O₂ Concentration channel.
12. Mean RER, which is the value for Mean on the RER channel. This is calculated by dividing the expired CO₂ value by the expired O₂ value to get the RER ratio as setup previously.
13. Record the T2-T1 value that is used to measure the mean values of CO₂, O₂ and RER.
14. Record the values in the Journal using one of the techniques described in Steps 6 or 7.
15. Record the values for the Mean CO₂ and O₂ concentrations in expired air in Table 1. Also record the RER and T2-T1 values.

Exercise 2: Changes in CO₂ and O₂ in a closed chamber, and RER in an ectothermic animal.

Aim: To determine the changes in CO₂ and O₂ volumes of an ectotherm and make a comparison with the values obtained from an endotherm.

Approximate Time: 15-20 minutes, depending on the animal.

Procedure

1. Use the same procedures used in Exercise 1 to record the CO₂, O₂ and RER values from an ectotherm.
2. Place the ectotherm (lizard, frog or snake) into the clean, dry small animal chamber and close the lid securely.
3. Mark the recording with comments that indicate the organism in the chamber.
4. Click Record to begin the recording.
5. Record until a concentration of 3% carbon dioxide is reached.

Note: Remove the animal from the chamber immediately in the CO₂ concentration reaches 3.5%.

6. Click Stop to halt the recording and Save your data file.
7. Remove the animal from the chamber and place it back in its container.

RER Setup and Data Analysis

1. Use the same procedures used in Exercise 1 to set up the RER channel showing a 1% change in CO₂ value.
2. Make any necessary changes on the RER channel by opening the RER Calculation dialog window.
3. Follow the procedures in Exercise 1 for data analysis of the beginning section of the recording for the ectotherm showing the step-like histogram.
4. Determine the oxygen (O₂), carbon dioxide (CO₂), and respiratory exchange ratio (RER) values. Note the T2-T1 value for the section of data with the step-like histogram.
5. Record the values for the Mean CO₂, O₂ concentrations, RER, and the T2-T1 value in Table 1.

Questions

1. During which experimental period was the endotherm's expired CO₂ the highest? In which period was it the lowest?
2. During which period was the endotherm's expired O₂ the highest? In which period was it the lowest?
3. During which period did the endotherm have the highest RER? In which period was the RER the lowest?
4. How do the values obtained from the endothermic animal compare to those from the ectothermic animal?
5. How do the RER values from these organisms relate to their overall metabolic rate?
6. Does ambient temperature have anything to do with their overall metabolism?
7. What would happen to the RER values for the ectotherm if you lowered the ambient temperature? The ectotherm?
8. For what reason should the animal not have eaten within an hour of performing these experiments? Evaluate the diet of these animals. How does diet correlate to the RER values?




Figure AMe-1-L4: The expired oxygen and carbon dioxide concentration and RER of a small animal as displayed in the Analysis window. The cursors are in position to measure the mean values of the selected interval.

Optional Exercise

RER is a temperature dependent value. As an optional exercise, the animals can be cooled and RER can be calculated as they return to normal body temperature or to room temperature.

Table AMe-1-L13: Table for recording Mean CO₂, O₂ and RER values for the endothermic and ectothermic organisms. T2-T1 values are also recorded.

Environmental Conditions		Organism	Mean Concentration (%) in the Chamber		Mean RER
Temperature (°C)	T2-T1 (min/sec) for RER calculation		CO ₂	O ₂	
		Endotherm			
		Ectotherm			

Weight Normalization Calculation

Normalized Expired O₂

$$\text{O}_2 / [(\text{weight(g)} / \text{mass unit})]^{\text{effective mass}}$$

Example: O₂ = 0.83 ml/min, weight = 25 grams,
mass unit = KG, effective mass factor (slope) = 0.75

$$\text{O}_{2\text{norm}} = 0.83 / (25/1000)^{0.75} = 13.2 \text{ ml/KG/min}$$

Using an effective mass factor of 1 eliminates any effective mass correction.

This can be repeated for the CO₂ values. When RER is calculated using normalized values, the weight of the organism is taken into account.