

8.12: Respiration Analysis

Introduction

The **Respiration Analysis Module** calculates physiologically relevant parameters from previously recorded respiration pressure and volume data.

This document includes a step by step tutorial for using the features of the *LabScribe Respiration Analysis Module*. To use the step by step guide, you will need a recording of recorded Respiration from any species.

Respiration Analysis: Step by Step


Offline Calculations

Sophisticated analysis can be done on a previously recorded Respiration data file using the **Offline Calculations** function of the **Respiration Advanced Analysis Module**. This analysis is performed using the offline **Respiration Calculations** dialog, so you should first become familiar with this dialog.

The Offline Respiration Calculations Dialog

To display the **Respiration Calculations** dialog and familiarize yourself with its features:

- 1) If it is not already open, open a respiration recording.
- 2) Select **Respiration Analysis** from the **Advanced** menu. This will open the offline **Respiration Calculations** dialog.
- 3) Familiarize yourself with the offline **Respiration Calculations** dialog, pictured below.
 - Across the top of the dialog, in the channel display area, you will see a sample of the raw data channel to be analyzed, including the selection between the two cursors. By default **Channel 1** is displayed.
 - On the left of the middle row are the tabbed dialogs used to configure the analysis.
 - At the right is the XY graph window in which the **Respiration Graph** is displayed.
 - Between the configuration dialogs and the graph are the editable lists of the **Cycles** to be analyzed and displayed.
 - Across the lower part of the dialog is the **Data Table** with the calculated average values for each of the analyzed groups of beats.




Respiration Calculations dialog.

To configure the analysis, the tabbed configuration panels at the left side of the middle row of the dialog are used.


To configure **Channels**:

- 1) Click the leftmost tab of the configuration dialogs, the one labeled **Channels**.
- 2) Choose the Volume and the Pressure Channel.



To configure the **Settings**:

- 1) Click the **Settings** tab, which is the second tab from the left in the configuration dialogs. The **Respiration Settings** configuration dialog will open.




- 2) Choose when the end of the breath should be determined.
- 3) Choose the point at which the Resistance should be calculated.
- 4) From the **Cycle Detection Threshold Sensitivity** menu, choose **2**. It is important that the cycle detection is set to the correct sensitivity. Adjusting the **Cycle Detection Threshold Sensitivity** number to higher numbers will lower the threshold at which a cycle is detected. Start at a low value; you will be able to adjust this later if you discover that cycles are being missed in the analysis.
- 5) Click the **Calculate** button just above the **Data Table** to start the analysis. The **Respiration Graph** will appear in the graph window at the right, and the **Data Table** will be populated with values.

Important: After any configuration settings are changed, click **Calculate** again, to trigger the revised analysis.

Once the **Channels** and **Settings** configuration dialogs are completed, it is possible to view the **Respiration Graph** and start the analysis.

To display the **Respiration Graph**:

- 1) The Respiration Graph should be automatically displayed in the XY Graph window, showing the selected cycle from the group of cycles specified in the **Cycles in Group** list to its left.
- 2) Use the menu indicated by the arrow at the lower left of the graph to **Copy Graph** to the clipboard, **Set the Y axis scale**, or **AutoScale the Y axis**.
- 3) Look at the **Respiration Graph** and familiarize yourself with its features.
 - Respiration parameters are indicated by the vertical blue Marks on the graph. The parameters that are shown are determined by the **Display** configuration dialog.
 - The specific cycle shown in the graph corresponds to the checked cycle in the **Cycles Selected** list to the left of the graph. The parameters and calculations from this cycle appear in the **Data Table**.
- 4) Change the cycle displayed by selecting a different cycle in the **Cycles Selected** list. The cycles are listed in order of their appearance in the data file.



Copy Graph
Set Y axis Scale
AutoScale the Y axis

The **Results** dialog includes basic information about the selection being analyzed,

- 1) Click on the **Results** tab to open the **Results** configuration dialog.
- 2) The **Results** dialog includes basic information about the selection being analyzed.

Channels	Settings	Results	
Volume Channel:		STPD Volume	
Pressure Channel:		Pressure	
Selection Start:		2011-02-10 00:02:01.670	
Selection Duration:	6.29	sec	
Selection Stop:		2011-02-10 00:02:07.959	

Data Table


All the data for each cycle, as well as the averaged values for all the cycles, is included in the **Data Table**.

*	Sel	Time	Rate	VE	VT	PIF	PEF	Insp. T	Exp T	FIT	Apneic T	Compliance	Resistance	Conductance
units		sec	bpm	Liters	Liters	L/sec	L/sec	sec	sec	#	sec	Liters/mmHg	mmHg.s/Liters	Liters/mmHg.s
1	Yes	4.420	13	8.282	0.618	-0.299	0.460	1.840	2.640	0.411	8.910	-127.937	0.283	3.538
2	Yes	8.900	15	7.917	0.513	-0.293	0.455	1.590	2.300	0.409	12.800	-262.193	0.273	3.665
3	Yes	12.790	15	7.938	0.534	-0.308	0.431	1.640	2.400	0.406	16.840	-295.662	0.255	3.925
5														
6	#		3	3	3	3	3	3	3	3	3	3	3	3
7	Mean		15	8.046	0.555	-0.300	0.449	1.690	2.447	0.408	12.850	-228.597	0.270	3.709
8	SD		1	0.168	0.045	0.006	0.013	0.108	0.143	0.002	3.238	72.477	0.012	0.161
9	Max		15	8.282	0.618	-0.293	0.460	1.840	2.640	0.411	16.840	-127.937	0.283	3.925
10	Min		13	7.917	0.513	-0.308	0.431	1.590	2.300	0.406	8.910	-295.662	0.255	3.538

The Respiration Data Table.

To use the **Data Table** and export values to the **Journal**:

- 1) Click **Table Options** at the bottom of the dialog to see a list of all the cardiac action potential parameters that can be displayed in the **Data Table**. These parameters and calculations are all defined in the **Algorithms** dialog, and are summarized below.
- 2) Choose the options you wish to include in the analysis and display in the **Data Table**. Choose whether you wish to display the **Time from the Start of the Selection** or the **Time of Day** of the recording. Click **OK**.
- 3) Click the asterisk in the upper left corner of the **Data Table**. The **Autosize** option adjusts the size of the cells for optimal display. The **Copy Selection** option will copy any selected cells to the clipboard.
- 4) Click **Algorithms** to see the definitions of the parameters and calculations. The definitions are also included below.
- 5) To copy all the calculated data in the **Data Table** to the clipboard, click the **Copy** button, or click the **Export** button to export the data. The data are exported in a tab (*.txt) or comma (*.csv) separated text file, and the graph can be exported as a Portable Network Graphics (*.png) or JPEG (*.jpg) image.
- 6) To load the analysis configuration for the current analysis, click **Save Template** to name and save the settings. Clicking **Load Template** when the module is reopened will display the list of previously saved templates.
- 7) Click **OK** to save the current configuration. The next time the offline **Respiration Calculations** dialog is opened, it opens with these settings.



Offline Calculation Algorithms: The offline calculations include:

- **Rate** : Respiration Rate (bpm)
- **VE**: Minute Vol. (VE) is calculated as $VT \times \text{Rate}$
- **VT**: Tidal Vol. Volume of air displaced between normal inspiration and expiration
- **PIF**: Peak Inspiratory Flow.
- **PEF**: Peak Expiratory Flow
- **Insp. T**: Inspiration Time
- **Exp. T**: Expiration Time
- **FIT**: Fractional Inspiratory Time, it is calculated as
$$\frac{\text{Inspiration Time}}{\text{Inspiration Time} + \text{Expiration Time}}$$
- **Apneic T**: Apneic Time. Time from end of breath to start of next breath
- **Compliance**: Dynamic Compliance
$$\frac{\Delta V}{\Delta P} \text{ at zero flow pt}$$
- **Resistance**: Dynamic Resistance
$$\frac{\Delta P}{\Delta F} \text{ at isovolume pts}$$
- **Conductance**: Dynamic Conductance
$$\frac{\Delta F}{\Delta P} \text{ at isovolume pts}$$